


# DIRECTIONS FOR USE: Correct Medication Labeling System

Some Components Patented & Patent Pending  
Pat. # US 6,955,002 B2; D538,851; D542,415


**1** WRITE IN APPROPRIATE INFORMATION USING A PERMANENT MARKER

<b>Afrin</b> Strength: _____ Expiration Date / Time: _____	<b>Marcaine</b> Strength: _____ % Expiration Date / Time: _____
<b>Afrin</b> Strength: _____ Expiration Date / Time: _____	<b>Marcaine <math>\pm</math> Epinephrine</b> Strength: <u>5 % / 1:200,000</u> Expiration Date / Time: <u>2-24 12:00</u>
<b>Antibiotic <math>\pm</math> 0.9% Saline Irrigation</b> Name and Strength: _____ Expiration Date / Time: _____	<b>Papaverine</b> Strength: _____ Expiration Date / Time: _____
<b>Contrast</b> Strength: _____ Expiration Date / Time: _____	<b>Contrast</b> Strength: _____ Expiration Date / Time: _____
<b>Lidocaine</b> Strength: <u>2 %</u> Expiration Date / Time: <u>2-24 12:00</u>	<b>Lidocaine</b> Strength: <u>2 %</u> Expiration Date / Time: <u>2-24 12:00</u>

**2** ADHERE LABEL ONTO SYRINGE




**3** PEEL PAPER BACKING HALFWAY BACK. ADHERE LABEL TO FLAG SURFACE. USE PRE-PRINTED OR MARK ON LABELS.




**4** PEEL PAPER BACKING OFF TO EXPOSE ADHESIVE. COLORED SYRINGE LABELS WILL COME OFF WITH BACKING.

**5** ADHERE FLAG TO OUTSIDE OF CUP. OPTIONAL: REMOVE ADHESIVE BACKING FROM BOTTOM OF CUP AND THEN ADHERE TO A SURFACE.



**6** PEEL OFF SYRINGE LABEL. ADHERE CORRESPONDING COLORED LABEL TO SYRINGE BARREL.




**7** SET SPECIMEN STRIP AND MARKING PEN ON BACK TABLE.


1 SPECIMEN ZONE	2 SPECIMEN ZONE	3 SPECIMEN ZONE	4 SPECIMEN ZONE
LEFT RIGHT	LEFT RIGHT	LEFT RIGHT	LEFT RIGHT



**REMEMBER TO TAKE A**



**8** PLACE SPECIMEN ON SQUARE. WRITE SPECIMEN NAME DIRECTLY ON PAD. USING PERMANENT INK.



OPTIONAL: IDENTIFY LEFT OR RIGHT. USE NUMBERS OR LETTERS TO COUNT SPECIMENS.

**SANDEL**  
the original exclusive sponsor of  
**AORN**  
Patient Safety First program.

**SANDEL** LLC  
Leaders in Healthcare Safety

Sandel Medical Industries, LLC  
19736 Dearborn Street, Chatsworth, CA 91311  
1-818-534-2500 www.sandelmedical.com