



DIRECTIONS FOR USE:
Z-Cord Holder™

Reorder No. 1117 & 1117-NS

- 1** PEEL BACK
ADHESIVE AND
PLACE IN DESIRED
POSITION.



- 2** INSERT TUBES
AND WIRES
BETWEEN
SLOTS.



The Z-Cord Holder is a disposable product.
Do not re-use.

Safety means NEVER having to say you're sorry.®

In the healthcare industry, the color orange is a trademark of Sandel Medical Industries, LLC

**SANDEL** LLC

Leaders in Healthcare Safety

Sandel Medical Industries, LLC
19736 Dearborn Street, Chatsworth, CA 91311
1-818-534-2500 www.sandelmedical.com