


DIRECTIONS FOR USE: DriFloor™

Reorder No. 2419, 2440, 2449

1

PLACE DriFloor™
IN AREAS WHERE
FLUID
COLLECTS.
WHITE SIDE
FACES
UP.



2

DISPOSE IN THE APPROPRIATE CONTAINER
BY LIFTING THE MAT IN A MANNER THAT IT
DOES NOT LEAK.

Caution: For safe use do not place on wet floor or in standing water. Do not over saturate to produce standing water. Use additional pads to absorb excess fluids.

SANDEL^{LLC}
Leaders in Healthcare Safety

Sandel Medical Industries, LLC
19736 Dearborn Street, Chatsworth, CA 91311
1-818-534-2500 www.sandelmedical.com

(SINGLE USE ONLY)