



DIRECTIONS FOR USE: - SINGLE TIP SKIN MARKER -

Pat. # US D538,851

Reorder No. 1019 (Sterile)
Reorder No. 1019-NS (Non-Sterile)



- 1. Remove purple cap for skin marker before use.**
- 2. Replace cap after use.**



Important!

1. For single patient use only.
2. Dispose of after use.
3. Skin surface must be clean and dried before marking the skin.
4. Alcohol or alcohol-based solutions will remove marks.

SANDEL
the original exclusive sponsor of
AORN
Patient Safety First program.


SANDEL LLC
Leaders in Healthcare Safety

Sandel Medical Industries
19736 Dearborn St. Chatsworth, CA 91311
818-534-2500