


# DIRECTIONS FOR USE: - SINGLE TIP PERMANENT MARKER -

Pat. # US D538,851

Reorder No. 1000-L (Sterile)  
Reorder No. 1000-N (Non-Sterile)



- 1. Remove black cap for permanent marker before use.**
- 2. Replace cap after use.**



Permanent Ink.

## Important!

1. Permanent Ink is not intended to be used on the skin.
2. For single patient use only.
3. Dispose of after use.

**SANDEL**  
the original exclusive sponsor of  
**AORN**  
Patient Safety First program.

**SANDEL**<sup>LLC</sup>  
Leaders in Healthcare Safety

Sandel Medical Industries  
19736 Dearborn St. Chatsworth, CA 91311  
818-534-2500