

DIRECTIONS FOR USE:
TIME OUT BEACON™

Reorder No. 811-NS

**Take Time
For
TIME OUT®**


Use over the instrument tray, back table or near the surgical site before incision.

The TIME OUT Beacon can be sterilized as follows, according to your protocol:

- Gamma
- ETO (Ethylene Oxide)
- Steam

The TIME OUT Beacon is a disposable product.

Do Not re-use.



SANDEL
the original exclusive sponsor of
AORN
Patient Safety First program.



Sandel Medical Industries, LLC
19736 Dearborn Street, Chatsworth, CA 91311
1-818-534-2500 www.sandelmedical.com